

# Minor Procedure Screen

## 1. Select the MP (Minor Procedure)

**MINOR OFFICE PROCEDURES**

Jessica A Medflow SEX/F AGE 30 BIRTH DATE 12/27/1983

Exam Date: 7/11/2014 Reason for:

Procedure:   require doctor to select ICD9

Exam Date	Reason	Procedure
7/30/2013		Botox Injection (add J0585), Bilateral Lower Lids, Bilateral Lower Lids
6/24/2013	Growth in and/or around ...	EYE OSISEE OP NOTE    Botox Injection (add J0585), Left Lower Lid
6/18/2013	Excessive Tearing, Exces...	Plugs Perm X, Left Lower Lid    Punctal Plugs 6 Month Dissolvable, Both Eyes
5/6/2013	Pain in and around eye	SEE OP NOTE
8/10/2012		Injection Sub Conj

**ICD9**

CORNEA/CONJUNCTIVA	LACRIMAL	EYELID	MISC
Foreign Body Removal - Superficial Foreign Body Removal - Embedded Foreign Body Removal - Cornea Subconjunctival Injection Keratotomy Incision Conjunctiva Drainage of Cyst Scraping Crn Dx F/Smr&/Culture Rmvl Crnl Epithe +-Chemocaut Corneal Relaxing Incision - Surgically Induced Ast Prokers Excision of Conjunctival Lesion	Punctal Plugs Temp 6 Month Dissolvable Punctal Plugs Temp 7 Day Dissolvable Punctal Plugs Permanent Dilation of Lacrimal Punctum w/ Irrigation Closure of Lacrimal Punctum - thermo	Epiletion Excision of Chalazion - Single Excision of Chalazion - Multiple Same Lid Excision of Chalazion - Multiple Diff Lid Foreign Body Removal from Eyelid Excision of Eyelid Lesion Closure of Eyelid by Suture Tarsorrhaphy Excision Benign Lesion or Margin (No Tag/Lid) Biopsy - Eyelid Intralesional Injection (up to 7)	Botox Injection (add J0585) Repair Retinal Detachment, Cryotherapy Revision of operative wound Paracenteses release aqueous Peripheral Iridotomy

- From the Minor Procedure screen, select the reason for the procedure by clicking in the box (marked REASON FOR) to bring up the pick list.

MINOR OFFICE PROCEDURES

Melinda K Peace SEX/F AGE 59 BIRTH DATE 12/22/1950

Exam Date 12/1/2010 Reason for COMBINED EXAM UPDATE SINGLE SCREEN

Procedure:  require doctor to select ICD9

Exam Date	Reason	Procedure
8/26/2010	Pain in and around eye	EYE OD/Removal Of FB, Conjunctiva, Slit Lamp, Right Upper Lid

- Choose an item off the pick list or free-type in the box under **Replacement Reason for Procedure** then hit **SELECT**

REASON FOR PROCEDURE Copyright Medflow, Inc. 1999-2010

Replace Current

Update From New

Append to Work

Undo Last Append

Clear Work

Clear All

Cancel Update

Replacement REASON FOR PROCEDURE

COMBINED EXAM UPDATE SINGLE SCREEN

- Double click on the procedure you are performing:

ICD9	DESCRIPTION
EYELID	Removal Of FB, Conjunctiva, Simple Removal Of FB, Conjunctiva, W/wo Slit Lamp Removal Of Foreign Body With Slit Lamp Biopsy Of Conjunctiva
LACRIMAL	Dilation Of Lacrimal Punctum, With Or Without Irrig: Punctum Dilation / Irrigation Removal Of Foreign Body, External Eye, Corneal, Closure Of The Lacrimal Punctum; By Thermocaut Punctum Plugs, Collagen Or Silicon Snip Incision Of Lacrimal Punctum Biopsy Of Lacrimal Gland
EYELID 11440 (078.10.216.1.)	Excision Benign Lesion Up To 0.5 Cm Excision Benign Lesion, 0.6 - 1.0 Cm Excision, Malignant Lesion, Lesion Diameter 0.5 or Excision, Malignant Lesion Diameter 0.6 - 1.0 Cm Biopsy Of Eyelid Epilation - Cautery Epilation - Forceps Excision Of Chalazion; Single Excision Of Chalazion; Multiple, Same Lid
EDM	Strabismus Surgery, Recession Or Resection Pro

- Choose the location of the procedure and hit **SELECT**

MINOR PROCEDURE Copyright Medflow, Inc. 1999-2010

Replace Current

Update From New

Append to Work

Undo Last Append

Clear Work

Clear All

Cancel Update

Replacement MINOR PROCEDURE

Excision Benign Lesion Up To 0.5 Cm


Current MINOR PROCEDURE

Excision Benign Lesion Up To 0.5 Cm

Push Up

5. Build the Op Note:

- a. Select the diagnosis from the upper left corner
- b. Select the Reason for Procedure from the bottom of the screen
- c. Select the Note associated with the procedure.



- d. Select the appropriate eye
- e. Select **Build Op Note**.

6. The completed op note will appear. You can free type to make changes/add/delete information. Once complete, **Put on Exam Document Only** or **Create a Separate Document**. (If you choose create separate document, it will bring up the document, click **Save/Exit** and it will go to the doctor's sign chart box.)

